


GRAPHICAL ABSTRACT



INTRODUCTION

As clinical proteomics applications start to reach maturity, this necessitates the requirement for robust and routine high throughput sample preparation workflows that allow processing of large sample cohorts. However, efficient sample preparation remains the *Achilles Heel* for mass spectrometry analysis, with current methods lacking the throughput, transferability and reproducibility required to deal with these large clinical sample numbers in a routinized laboratory setting. To address these we focus on the implementation of versatile and automatable magnetic bead based sample preparation. Magnetic beads are considered desirable since these are easy to handle, linearly scalable, and high throughput compatible with the relatively simple integration of a magnetic stand in a variety of liquid handling stations making it fairly agnostic for the liquid handling or magnetic bead handling station that may be present in the sample preparation laboratory. We have previously demonstrated protein and peptide clean-up workflows using magnetic HILIC for Solid Phase Extraction (HILIC SPE) from a broad range of common contaminants, while the use of magnetic beads for highly efficient phosphopeptide enrichment has been demonstrated (Tape *et al.* 2014, Baath *et al.* 2018). In the current work we demonstrate the coupling of these steps into an automated workflow for phosphoproteome profiling requiring minimal liquid handling.



Although our aim is to fully automate mass spectrometry workflows, the protocols are also suitable for manual preparation with the ability to perform parallel sample processing using a magnetic stand. The protocols can be transferred to a variety of liquid or bead handling systems. All experiments were performed on a KingFisher™ Duo (ThermoFisher) magnetic bead handling station (protocols available on request).




LC MS/MS & BIOINFORMATICS

Samples were analysed using an AB SCIEX TripleTOF 6600 coupled to a Dionex nanoRSLC via a nanoSpray III interface with 60 minute gradient. Spectral data was searched using PEAKS Studio 6 (Ma *et al.*, 2003, Bioinformatics Solutions Inc). In the case of HCT 116 samples a Swiss-Prot mammalian database, supplemented with sequences of common contaminant proteins. A 0.1% and 1% FDR cut-off was applied at the PSM and peptide/protein levels respectively.

OPTIMIZATION BLOCK 1: METHODS, RESULTS & DISCUSSION


PROTEIN CLEAN UP, REDUCTION, ALKYLATION AND DIGESTION

AIM: To evaluate particle chemistry, and the effect of on-bead reduction and alkylation on automated sample processing.



REFERENCES


- Bath ST, Tollenaere MAX, Ruther PL, Gonzalez-Franquesa A, Prabhakar BS, Bekker-Jensen SH, Deshmukh AT, Olsen JV. 2018. Protein aggregation capture on microparticles enables multi-purpose proteomics sample preparation. *bioRxiv* 447904; doi: <https://doi.org/10.1101/447904>
- Ma B, Zhang K, Hendrie C, Liang C, Li M, Doherty-Kirby A & Lajoie G. 2003. PEAKS: powerful software for peptide de novo sequencing by tandem mass spectrometry. *Rapid Comm Mass Spec.* 17(20), 2337-2342.
- MacLean B, Tomazela DM, Shulman N, Chambers M, Finney GL, Frewen B, Kern R, Tabb DL, Liebler DC, MacCoss MJ. 2010. Skyline: an open source document editor for creating and analyzing targeted proteomics experiments. *Bioinformatics.* 26(7):966-8.
- Tape CJ, Worboys JD, Sinclair J, Gourlay R, Vogt J, McMahon KM, Trost M, Lauffenburger DA, Lamont DJ, Jørgensen C. 2014. Reproducible automated phosphopeptide enrichment using magnetic TiO₂ and Ti-IMAC. *Analytical Chemistry.* 86 (20), 10296-302.
- Thingholm TE, Larsen MR. 2016. Phosphopeptide Enrichment by Immobilized Metal Affinity Chromatography. *Methods. Mol Biol.* 1355:123-33.




OPTIMIZATION BLOCK 3 & INTEGRATION WITH BLOCK 1: METHODS, RESULTS & DISCUSSION

PHOSPHOPEPTIDE ENRICHMENT


AIM: Evaluate compatibility of phosphoenrichment strategies with clean-up of sample prior to MS analysis




In order to optimize Block 3 we evaluated a range of phosphopeptide enrichment tools with particular interest in enrichment efficiency, capacity and recovery. To this end a range of functionalities were assessed including magnetic beads having immobilized titanium ions (Ti-IMAC), zirconium ions (Zr-IMAC), and titanium and zirconium oxides (TiO₂ and ZrO₂), using a standard protein mix. All enrichment protocols were used as per manufacturers instructions (ReSyn Biosciences).




A Casein & BSA tryptic digest was processed (automated on a KingFisher™ Duo), using the range of magnetic micro-particles. Post phosphopeptide enrichment samples were spiked with 3x isotopically labelled peptide standards correlating to the phosphopeptides of casein, VPQLIVPN(pSer)AEER, YKVPQLIVPN(pSer)AEER and FQ(pSer)EEQQQTDEQDK. MALDI-TOF MS1 spectra were generated and the ratio of light to heavy (L:H) phosphopeptide was calculated for each sample to determine enrichment efficiency. Experiments were performed in triplicate with each experiment further analysed in triplicate. The Ti-IMAC and Zr-IMAC enriched samples showed similar L:H ratios followed by TiO₂ and last ZrO₂. In the case of ZrO₂, 10x more beads had to be applied as compared to the other 3 enrichment chemistries. It is apparent that the TiO₂ may further require a higher quantity for improved enrichment. Considering the potential complementarity of the chemistries by Tape *et al.* in 2014. Future work will include mixing the range of functionalities with the aim of increasing phosphoproteome coverage.



Samples of 200 µg protein extracts in 1% SDS were processed using FASP (TEAB as digest buffer) or HILIC (TRIS, AmmBic, PBS or TEAB as digest buffers) using a magnetic handling station. Post digestion 5µl from each sample was kept for LC-MSMS analysis while remainder was used for phosphopeptide enrichment on Ti-IMAC magnetic micro-particles. There was a 10% drop in peptide and protein identifications in the case of the FASP processed samples and 15% drop-off when on-bead digestion was performed in PBS with the HILIC-based protocol. The peptide recoveries are relative to 10977 and proteins ID's to 1947 (Peptides: HILIC SPE in TRIS; Proteins HILIC SPE in AmmBic). Experiments were performed in duplicate.



The specificity of the Ti-IMAC enrichment (A) appeared unaffected by the selection of tryptic digest buffer nor the method of digest preparation (i.e. FASP vs HILIC SPE), with 99% of the identified peptides being phosphorylated. The selectivity (B) for type of phosphorylation was further not significantly affected by buffer or method of preparation. As expected approximately 80% of the identified phosphopeptides were pSer, followed by 20% pThr and 2-4% pTyr.




The distribution of phosphorylation levels i.e. single vs double vs triple phosphorylation, was significantly affected by the buffer used for digestion (D, E, F & G), and to a lesser extent the method of digest preparation (i.e. C vs D FASP vs HILIC SPE). The combination of HILIC SPE with digestion in AmmBic (E) resulted in predominantly singly phosphorylated peptides being identified. On-bead digestion performed in the presence of PBS yielded a significantly higher amount or multiply phosphorylated (double and triple) peptides, perhaps stabilised by the addition of native phosphate buffer, while the FASP-TEAB (C), HILIC-TEAB (D) and HILIC-TRIS (G) showed a more even distribution between singly and doubly phosphorylated peptides. This phenomenon will be assessed in future experiments where the HILIC based workflow is expanded for applications using multiple or complementary buffers e.g. on-bead digestion in TRIS followed by PBS incubation of HILIC micro-particles prior to phosphopeptide enrichment.

ALTERNATE BLOCK 3: METHODS, RESULTS & DISCUSSION


SAX CLEAN-UP PRE-MS ANALYSIS

AIM: Evaluate the application of SAX for clean-up of peptides for direct MS analysis of peptides



With the aim of fully automating the sample preparation workflow, we evaluated the option of using strong anion exchange (SAX) as a clean-up step prior to MS analysis. The use of this chemistry has several advantages for clean-up, primarily that binding occurs at alkaline pH (i.e. compatible with elution after Ti-IMAC), and elutes at acidic pH thereby removing the necessity for drying of samples (required with HILIC SPE of peptides).

A tryptic peptide digest of 10 µg was processed using magnetic HILIC SPE or SAX SPE using automated protocols. SAX SPE was performed with binding in 50 mM ammonium formate, and elution in 0.2% formic acid. HILIC SPE for peptides was performed as described previously. Experiments were performed in duplicate. Initial evaluation revealed that peptide and protein identifications were significantly reduced in the case of SAX processed samples (A) indicating loss of sample under the binding and elution conditions tested. The peptide recoveries are relative to 10393 peptide and 2177 protein ID's (HILIC SPE). The peptide GRAVY (B) and pI distributions (C) indicated that hydrophobic peptides and peptides with alkaline pI values did not bind well during SAX SPE, providing an explanation for the reduced binding efficiency. Due to the potential to go directly to MS without offline drying, and the fractionation potential for deep proteome profiling, we intend to elaborate on the current SAX workflow, and evaluate a range of binding and elution conditions.



CONCLUSIONS & FUTURE WORK

- This work describes a pipeline for phosphoproteome profiling that allows for automated processing of up to 96 samples in parallel in less than 8hrs with all methods automated on a magnetic bead handling station, KingFisher™ (protocols available on request)
- The pipeline consists of modular and interchangeable BLOCKS including:
 - Efficient protein isolation in detergents or denaturants such as SDS or UREA with integrated reduction and alkylation, removal of contaminants, and on-bead tryptic digestion of proteins on magnetic HILIC or SAX microparticles
 - Phosphopeptide enrichment using magnetic Ti-IMAC, Zr-IMAC, TiO₂ or ZrO₂ that does not require sample desalting or lyophilization post digestion allowing for seamless and automatable protocol integration
 - Optional desalting of peptides using HILIC or SAX magnetic micro-particles (with Ti-IMAC or Zr-IMAC this may be omitted)
- Future work will focus on:
 - Improving the on-bead reduction and alkylation protocol
 - Optimization of SAX SPE for protein recovery and clean-up with coupling to on-bead digestion and peptide fractionation for deeper proteome and phosphoproteome profiling
 - Evaluate multiple buffer elution from HILIC SPE with coupling to the phosphopeptide enrichment protocol
 - Combination of various phosphopeptide enrichment products to improve phosphoproteome coverage (eg. Ti-IMAC and Zr-IMAC)
 - Optimise SAX SPE for peptide binding, desalting, and concentration of phosphopeptides directly from alkaline elution, allowing for direct MS analysis of samples