



our future through science

Development of a fully automated magnetic workflow for phosphoproteome profiling

Stoyan Stoychev^a, Previn Naicker^a, Iresyhn Govender^a, Siphso Mamputha^a, Sindisiwe Buthelezi^a, Isak Gerber^a, Justin Jordaan^b
^a Council for Scientific and Industrial Research, Pretoria, South Africa
^b Rhodes University, Grahamstown, South Africa
E-mail: ssoychev@csir.co.za



Poster QR Code

GRAPHICAL ABSTRACT

OVERALL AIMS



- HT COMPATIBLE
- AUTOMATED
- VERSATILE WORKFLOWS
- REAGENT COMPATIBILITY
- VENDOR INDEPENDENT
- MODULAR

CHALLENGES

- AUTOMATION
- INTEGRATION
- REPRODUCIBILITY
- DATA QUALITY

ADVANTAGES

- HIGHLY REPRODUCIBLE
- LINEARLY SCALABLE
- 96 SAMPLES in <8 HRS
- 5 MIN PER SAMPLE



GENERAL PROTEOMICS WORKFLOWS

PROPOSED MAGNETIC PROTEOMICS WORKFLOWS

Current Work:

- Optimization of HILIC and SAX for automated protein clean-up
- On and off-bead reduction and alkylation

Future Work:

- Use of LysC & Trypsin for on-bead digestion.
- In-depth evaluation of SAX for sample clean-up

Future Work:

- Evaluate tools for fractionation including HILIC, SAX, Hydrizide and IP
- Integrate and automate fractionation strategies
- Evaluate tools for proteome profiling

Previous Work:

- Automation of phosphopeptide enrichment
- Evaluation of HILIC for protein and peptide clean-up

Current Work:

- Integration of HILIC & Ti-IMAC workflow under range of digest conditions & comparison to FASP & Ti-IMAC workflow

Future Work:

- Evaluation of workflows for complex samples
- Optimization of SAX for peptide desalting

Aims of Current Work:



- Optimization and automation of sample clean-up prior to phospho-enrichment
- Integration and automation of sample clean-up and phosphopeptide enrichment
- Evaluation of available tools for improving global phosphoproteome coverage, preliminary investigation on compatibility of enrichment chemistry

INTRODUCTION

As clinical proteomics applications start to reach maturity, this necessitates the requirement for robust and routine high throughput sample preparation workflows that allow processing of large sample cohorts. However, efficient sample preparation remains the *Achilles heel* for mass spectrometry analysis, with current methods lacking the throughput, transferability and reproducibility required to deal with these large clinical sample numbers in a routine laboratory setting. To address these we focus on the implementation of versatile and automatable magnetic bead based sample preparation workflows. Magnetic beads are considered desirable since these are easy to handle, linearly scalable, and high throughput compatible with the relatively simple integration of a magnetic stand in a variety of liquid handling stations making it independent of the liquid handling or magnetic bead handling that may be present in a sample preparation laboratory. We have previously demonstrated protein and peptide clean-up workflows using magnetic HILIC for Solid Phase Extraction (HILIC SPE) from a broad range of common contaminants. Further, the use of magnetic Ti-IMAC for highly efficient phosphopeptide enrichment has previously been illustrated by Tape *et al.* 2014 and Baath *et al.* 2019. In this study we demonstrate the coupling of the clean-up to phosphopeptide enrichment. We further evaluate the option to combine a range of phosphopeptide enrichment chemistries (Ti-IMAC, Zr-IMAC, TiO₂, ZrO₂, Fe-NTA, Fe-IMAC) for possible deeper phosphoproteome coverage, by identifying buffers suitable for enrichment using various combinations of bead chemistries.

AUTOMATION

Although our aim is to fully automate mass spectrometry workflows, the protocols are also suitable for manual preparation with the ability to perform parallel sample processing using a laboratory magnetic stand. The protocols can be transferred to a variety of liquid or bead handling systems. All current experiments were automated on a KingFisher™ Duo (ThermoFisher) magnetic bead handling station (protocols available on request).




Automation of magnetic clean-up workflows can be achieved with KingFisher magnetic bead handling stations. Examples of the conditions used for HILIC SPE of proteins using ACN (left) with SAX clean-up using alkaline conditions (middle), and phosphopeptide enrichment (right). For the clean-up workflows, the on-bead digestion takes place at position 8 since this is the position of the heating element. Beads may be transferred to position 2 or 3 if an extra elution step is required, or on-bead reduction and alkylation are included. Similarly, an extra elution step may be included for phosphopeptide elution if required.

LC MS/MS & BIOINFORMATICS



Samples were analysed using an AB SCIEX TripleTOF 6600 coupled to a Dionex nanospray II interface with 60 minute gradient. Spectral data was searched using PEAKS Studio 5 (Ma *et al.*, 2003, Bioinformatics Solutions Inc). In the case of HCT 116 samples a Swiss-Prot mammalian database, supplemented with sequences of common contaminant proteins. A 0.1% and 1% FDR cut-off was applied at the PSM and peptide/protein levels respectively.

OPTIMIZATION OF SAMPLE CLEAN-UP

AIM: To evaluate particle chemistry, and the effects of on-bead reduction and alkylation on automated sample processing prior to phosphopeptide enrichment.



Magnetic beads with HILIC or Strong Anion Exchange (SAX) chemistry were evaluated for the automation of protein clean-up, with on-bead reduction and alkylation (B) (off-bead run as a control (A)), followed by on-bead trypsin digestion. We have previously optimized digestion to 4 hrs using Trypsin (1:10 in 20 mM AmmBic pH 8). Due to the potential compatibility to extend the range of contaminants that can be removed by automated clean-up, we evaluated the use of SAX (with Urea as an example contaminant) as an alternative for clean-up of proteins using HILIC SPE. Proteins were extracted from HCT 116 cells using 8M Urea, and 10 µg aliquots were isolated using magnetic HILIC SPE or SAX SPE. Briefly, for HILIC SPE, proteins were adsorbed to magnetic HILIC microparticles using acetonitrile (ACN - 15%) under acidic conditions (ammonium acetate, pH 4.5), and washed with 95% ACN. For SAX SPE, proteins were adsorbed to SAX in Formate or Borate buffer (50mM) at alkaline pH of 10, and washed twice with the same buffer. All experiments were performed in duplicate.





SUMMARY

- HILIC and SAX are suitable for automated protein isolation and on-bead digestion
- Fully automated on-bead reduction and alkylation can be performed, but some optimization of this step is still required to improve sample coverage
- SAX shows promise for sample clean-up, with the possibility to extend the range of contaminants that can be removed in the fully automated workflow.

INTEGRATION OF HILIC CLEAN-UP WITH PHOSPHO-ENRICHMENT: Ti-IMAC

AIM: Evaluate compatibility of clean-up with the phosphoenrichment prior to MS analysis






SUMMARY

- The automated clean-up and digestion of proteins was successfully coupled to phosphopeptide enrichment without the requirement for additional desalting and/or concentration steps
- Automated HILIC SPE with on-bead trypsin digestion performed favourably in comparison to FASP with in-solution protein digestion
- The selection of digestion buffer significantly affects the identification of singly vs multiply phosphorylated peptides, this warrants further evaluation to potentially improve phosphoproteome coverage


OPTIMIZATION OF PHOSPHOPEPTIDE ENRICHMENT BUFFERS

AIM: Identify the ideal buffer conditions for phosphopeptide enrichment for a range of enrichment chemistries, for future integration into automated workflow for increased phosphoproteome coverage.

We evaluated a range of phosphopeptide enrichment tools with particular interest on how buffer composition affects phosphopeptide enrichment efficiency and recovery. To this end we assessed a range of commercially available magnetic bead chemistries (ReSyn Biosciences), including Ti-IMAC (chelated Ti³⁺ ions), Zr-IMAC (Zr⁴⁺ chelated ions), TiO₂ and ZrO₂ (titanium and zirconium dioxide nanoparticles), and two prototype magnetic beads containing Fe³⁺ ions chelated to two different supports (Fe-IMAC and Fe-NTA, ReSyn Biosciences). For initial evaluation of buffer composition a relatively simple mixture was generated from a tryptic digest of Casein and BSA. The optimal conditions will be applied to complex lysates (work currently in progress). For all experiments 1 mg of beads was used for enrichment, except in the case of ZrO₂ where 10 mg was evaluated (high density nanoparticles attached to the polymer beads).




Due to possible increased competition in complex samples, these conditions will be investigated further. Unlike Ti-IMAC, for Zr-IMAC high ACN reduced specificity (BC7). For TiO₂ a low concentration of GA (0.1M, BC5) improved recovery with no loss in specificity, while an increased concentration (0.5M, BC6) reduced recovery. The increase in ACN concentration (BC7) improved recovery but reduced specificity. For ZrO₂, the addition of 0.1M GA improved both specificity and recovery (BC8). Similar to TiO₂, increased ACN (BC10) improved recovery but reduced specificity.



BC	Binding buffer	Wash buffer 1	Wash buffer 2
5	50%ACN/0.1% Acetic, 0.1M GA	80% ACN, 1% TFA	10% ACN, 0.2% TFA
11	80%ACN/0.1% Acetic, 0.1M GA	80% ACN, 1% TFA	10% ACN, 0.2% TFA
12	80% ACN, 5% TFA, 0.1M GA	80% ACN, 1% TFA	10% ACN, 0.2% TFA
13	80% ACN, 5% TFA, 0.1M LA	80% ACN, 1% TFA	10% ACN, 0.2% TFA
14	80%ACN/0.1% Acetic, 0.1M GA	80% ACN, 1% TFA	10% ACN, 0.2% TFA

The aim of experiment 3 was to further evaluate the effects of hydroxy acids including glycolic (GA), lactic (LA) and tartaric (TA) acids on phosphopeptide enrichment with Ti-IMAC, Zr-IMAC, TiO₂ and ZrO₂ magnetic microparticles (Sugiyama *et al.*, 2007). In all cases (except ZrO₂) the addition of TA (BC13 and BC14) resulted in significant decrease in phosphopeptide recovery. The decrease of GA from 1 to 0.1M GA improved phosphopeptide recovery without major effect on specificity (refer BC1 vs BC11 for each chemistry) while 0.1M LA performed similarly to 0.1M GA (BC12 vs BC11). In the case of Ti-IMAC and Zr-IMAC the lower pH buffer BC11 (TFA with 0.1M GA) performed better than BC5 (Acetic with 0.1M GA). This was opposite when comparing TiO₂ and ZrO₂, with TiO₂ showing an increase in recovery (BC5 vs BC11). The combined results from three experiments were used to derive a buffer table to guide future work where combinations will be applied to complex lysates.



Buffer Compatibility Table		
Enrichment Chemistry	Preferred Buffers	
Ti-IMAC	BC1	BC12
Zr-IMAC	BC11	BC12
TiO ₂	BC5	BC11
ZrO ₂	BC5	BC11
Fe-NTA	BC2	
Fe-IMAC	BC2	
MIXED *	BC5	BC11

* - Possibly suitable for mixture of chemistries

CONCLUSIONS & FUTURE WORK

- This work aims to develop a fully automated magnetic bead based workflow for deep phosphoproteome profiling that allows for processing of up to 96 samples in parallel in less than 8hrs
- The pipeline consists of modular and interchangeable blocks including:
 - Efficient protein isolation from detergents or denaturants (including SDS or UREA) with integrated reduction and alkylation, and on-bead tryptic digestion of proteins on magnetic HILIC or SAX microparticles
 - Phosphopeptide enrichment using magnetic Ti-IMAC which does not require sample desalting or lyophilization post digestion, enabling seamless and automatable protocol integration
- To improve global phosphoproteome coverage we have identified buffer conditions suitable for combining phosphopeptide enrichment chemistries. We are in process of applying these combinations to complex samples for identifying ideal conditions.
- The features of this enrichment include:
 - All chemistries are from the same supplier negating the possible effect of the support material.
 - Simple integration into automated protocols owing to their magnetic format.
- Future work will focus on:
 - Improving the on-bead reduction and alkylation protocol
 - Optimization and evaluation of alternate methods of clean-up such as SAX SPE and PAC (Baath *et al.*, 2019)
 - Evaluate a range of elution buffers for HILIC SPE for coupling to phosphopeptide enrichment.
 - Application of combinations of phosphopeptide enrichment chemistries for improved phosphoproteome coverage
 - Optimization of off-bead peptide fractionation (prior to enrichment) for deeper phosphoproteome coverage

REFERENCES

- Bath ST, Tollenaar MA, Ruther PL, Gonzalez-Franquesa A, Prabhakar BS, Bekker-Jensen SH, Deshmukh AT, Olsen JV, 2019. Protein aggregation capture on microparticles enables multi-purpose proteomics sample preparation. Mol Cell Proteomics. mcp.T1R118.001270
- Ma B, Zhang K, Hendrie C, Liang C, Li M, Doherty-Kirby A & Lajoie G. 2003. PEAKS: powerful software for peptide de novo sequencing by tandem mass spectrometry. Rapid Commun Mass Spec. 17(20), 2337-2342.
- Maclean B, Tomazela DM, Shulman N, Chambers M, Finney GL, Freuden B, Kern R, Tabb DL, Liebler DC, MacCoss MJ. 2010. Skyline: an open source document editor for creating and analyzing targeted proteomics experiments. Bioinformatics. 26(7):966-8.
- Sugiyama N, Masuda T, Shinoda K, Nakamura A, Tomita M, Ishihama Y. 2007. Phosphopeptide Enrichment by Aliphatic Hydroxy Acid-modified Metal Oxide Chromatography for Nano-LC-MS/MS in Proteomics Analysis. Mol Cell Proteomics. 6(1):313-9.
- Tape CJ, Watsboys JD, Sinclair J, Goussay R, Vogt J, McNaughton KM, Trost M, Lauffenburger BG, Lamont DJ, Jørgensen C. 2014. Reproducible automated phosphopeptide enrichment using magnetic TiO2 and Ti-IMAC. Analytical Chemistry. 86 (20), 10296-302.
- Thingholm TE, Larsen MR. 2016. Phosphopeptide Enrichment by Immobilized Metal Affinity Chromatography. Methods. Mol Biol. 1355:123-33.

