

AUTOMATED SAMPLE PREPARATION FOR BOTTOM-UP URINARY PROTEOME PROFILING IN CLINICAL PROTEOMICS

Ireshyn Govender^{a,b}, Previn Naicker^a, Andrea Ellero^c, Rethabile Mokoena^{a,c}, Stoyan Stoychev^b

^a Next Gen Health, CSIR, Pretoria, South Africa, ^b ReSyn Biosciences, Pretoria, South Africa, ^c University of Pretoria, Pretoria, South Africa, ^d University of the Witwatersrand, Johannesburg, E-mail: igovender@resynbio.com

GRAPHICAL ABSTRACT


COLLECTION OF CLINICAL URINE SAMPLES



NOVEL AUTOMATED URINE SAMPLE PROCESSING



DIFFERENTIATION IN CLINICAL COHORT









INTRODUCTION

- Urine is a diverse source of information related to a patient's health and is gradually becoming a go-to biofluid for clinical proteomics.
- There is currently no standard operating procedure for automated, reproducible & robust sample processing of urine for mass spec-based clinical proteomics
- Here, we present a novel workflow with direct protein capture on HILIC magnetic beads for Urine Low-volume Clean-up with on-bead digestion (H μ LC), and compare this with a published on-membrane (OM) method developed by Berger *et al* (2015), as amended by Winter *et al* in 2021.
- The primary advantages of H μ LC is that it only requires 100 μ l of urine, with direct on-bead processing and full automation.
- The H μ LC workflow was applied in two pilot studies: [A] Pancreatic Duct Adenocarcinoma (PDAC) & [B] Acute Kidney Injury (AKI)
- For Study A: Urine was collected from patients with resectable (RST) and irresectable or metastatic (iRST/M) PDAC at Chris Hani Baragwanath Hospital.
- Urine from AKI+/- patients was collected at the Tshepong Hospital Complex.
- We aimed to map the urinary proteomes of PDAC and AKI patients against a healthy cohort to gain clinical insights from their proteome profiles.


METHODS

- The OM (on-membrane) method was performed according to Winter *et al* (2021)
- The 2 workflows were compared using a 3 by 3 approach (n = 9 tech replicates). Peptides were analysed using an Ultimate 3000 coupled to a TripleTOF 5600.
- The H μ LC workflow was applied in two distinct pilot studies (AKI and PDAC) to evaluate variations in the protein profiles between the cohorts.
- Data was searched using a study specific library in Spectronaut 17 and differentially abundant proteins (DAPs) were analysed using GO annotation software.


HULC WORKFLOW: LIQUID HANDLING



Automated processing on KF Flex, protocol available from ReSyn Bio on request




Pierce Quantitative Colorimetric Peptide Assay



75 μ m ID, 25 cm, CSH column, nanoUPLC 30 min gradient at 300nL.min⁻¹, 48VW SWATH

Instrument, plate and sample images courtesy of biorender.com



SIMPLER


H μ LC: Direct on-bead urinary proteomics workflow: bind – wash – digest – LCMS

BETTER

Increased urinary proteome coverage with high reproducibility from low urine volume (100 μ l)


FASTER

Rapid magnetic workflow without time consuming precipitation or centrifugation, ~2.5 min per sample







SIMPLER BETTER FASTER

RESULTS: BENCHMARKING OF WORKFLOW



- A: The H μ LC method showed greater total peptide yield and lower CV's for peptide and protein groups than the OM method.
- B: There was a large overlap between the two methods as illustrated, with 7711 and 7477 peptides identified in the H μ LC and OM workflows respectively, corresponding to 1140 and 1069 protein identifications
- C: Minor differences were observed in MW, pI and GRAVY score for the unique proteins identified in each workflow
- D: PCA plot showed that samples prepared using OM had a greater spread of data, indicative of poorer reproducibility between the technical replicates, while results from samples prepared with H μ LC were more tightly clustered.
- E: The H μ LC method appeared to identify a higher proportion of lower abundant proteins (16% vs 12% for OM)



RESULTS: PILOT EVALUTION

PDAC PILOT

- A: A total of 78 DAPs were identified in iRST/M vs RST groupings, at 1% FDR, ≥ 2 fold change. The major biological pathways represented in this data were the immune system, complement cascade, and proteins associated with hemostasis
- B: Clustering analysis revealed some overlap between patients with RST and iRST/M.

AKI PILOT

- C: For AKI samples, 137 DAPs were identified with the same parameter settings. The proteins RBP4, CYTC and B2M were found in increased abundance with AKI+, and are known to be associated with kidney damage (Vaidya *et al.*, 2008)
- D: The AKI + and - cases formed distinct clusters during PCA analysis

CONCLUSIONS

- We have developed a novel workflow, H μ LC, suitable for the low-volume direct automated processing of clinical urine samples without the need for centrifugation or precipitation.
- Automated sample processing of 96 samples on KF Flex at <2.5 min per sample
- Clinical urine samples processed with H μ LC could successfully distinguish patients with AKI, and resectable vs irresectable/metastatic PDAC cases.
- The workflow may find clinical applications for the analysis of urine samples in high-throughput

REFERENCES

- Berger ST *et al.*, 2015. Blotting-High Throughput Polyvinylidene Fluoride (PVDF) Membrane-Based Proteomic Sample Preparation for 96-Well Plates. *Mol Cell Proteomics*. 14(10):2814-23. doi: 10.1074/mcp.O115.049650.
- Bernhardt OM *et al.*, 2012. Spectronaut A fast and efficient algorithm for MRM-like processing of data independent acquisition (SWATH-MS) data. *Proceedings of the 62th ASMS Conference on Mass Spectrometry and Allied Topics*, 2012, Vancouver, BC, Canada.
- Vaidya VS *et al.*, 2008. Biomarkers of acute kidney injury. *Annu Rev Pharmacol Toxicol*. 48:463-93. doi: 10.1146/annurev.pharmtox.48.113005.094615
- Winter S *et al.*, 2021. Urinary proteome profiling for stratifying patients with familial Parkinson's disease. *EMBO Mol Med* 13: e13257. doi: 10.15252/emmm.202003327

CONFLICT OF INTEREST DISCLOSURE

Stoyan Stoychev, Andrea Ellero and Ireshyn Govender are employed by ReSyn Biosciences, proprietors of MagReSyn[®] technology.